

**Patient**

Last name: \_\_\_\_\_

First Initial: \_\_\_\_\_

Patient ID: \_\_\_\_\_  N  W


Bilateral  Left only  Right only

Practitioner: \_\_\_\_\_ Please change to your name if necessary.

Facility: \_\_\_\_\_

**Construction • Features • Options**

**FLOOR TO HEIGHT MEASUREMENTS** All measurements must be in mm




**Shell Configuration**

- Solid Ankle
- Reinforced
- Posterior Stop
- Attach to Shoe
- Double Action
- Tamarack **Please Select Option**
- Neutral  75  85  95
- Dorsi Assist  75  85  95

**Construction • Features • Options**

All measurements must be in mm



**LEGEND**

- A1: Brace Length
- A2: Brace Height
- A3: Calf Circumference
- A4: Calf ML
- A5: Calf Height
- A6: Narrowest Ankle Circ.
- A7: Narrowest Ankle ML
- A8: Narrowest Ankle Height
- B1: Ankle ML
- B2: Ankle Height
- B3: 1st MET AP
- B4: 1st MET Circ.
- B5: 1st MET ML
- B6: Mid-Arch Circumference
- B7: Mid-Arch AP
- B8: Diagonal Heel Circ.
- B9: Diagonal Heel AP

**Cast Correction**

**SAGITTAL ANKLE ALIGNMENT** (Dorsiflexion-Plantarflexion)

- Correct to 3-4° DF
- Correct to \_\_\_\_\_°  DF  Do not correct
- PF  Cast alignment OK

**SAGITTAL KNEE CORRECTION**

- Correct to neutral
- Correct to \_\_\_\_\_°  Flexion  Do not correct
- Extension

**CORONAL HINDFOOT ALIGNMENT**


- Correct to vertical
- Correct to \_\_\_\_\_°  Varum  Do not correct
- Valgum

**CORONAL KNEE CORRECTION**

- Correct to neutral
- Correct to \_\_\_\_\_°  Varum  Do not correct
- Valgum

**FOREFOOT ALIGNMENT** NOTE: Drawings show finished orthosis.

Choose forefoot alignment. Write posting height - in mm - if needed.



Toe Out: Toe Out Angle: (+) To Out (-) Toe In

Left: \_\_\_\_\_

Right: \_\_\_\_\_

Ankle Joint Alignment  Mechanical  Other  Anatomical

Inner Boot:  Polyethylene  Softy Foam (white only)  Duraflex\*

**Standard**

- Add extra navicular padding (boney pronators only)
- Add plastizote to malleoli

Outer Frame:  Full-length under plantar surface  Trimmed distal to met. heads under plantar surface  Trimmed just proximal to met. heads under plantar



## Pads and Straps

Straps: **Standard** Add toe abduction strap  
(tibial & instep straps)  Add forefoot strap

Strap Color: **White is Standard** \_\_\_\_\_

Transfer Pattern:\*

Padding Color: **White is Standard** \_\_\_\_\_

### LINERS / PADDING / INSERTS

- None **Standard**
- Plantar Footplate
- Calf Section
- Thigh Section

## Special Instructions

Rush Order

## More Shell Configuration Items

### TONGUE

- Calf Section
- Thigh Section

### MATERIALS

- Polypropylene **Standard**
- Copolymer

### THICKNESS

- 1/8
- 3/16
- 5/32
- 1/4

### KNEE CONTROL

- None **Standard**
- Medial Condylar Extension
- Lateral Condylar Extension

## Component Selection

### OTHER ANKLE JOINT OPTIONS

Standard Action

- No motion
- Stop Motion \_\_\_ ° DF \_\_\_ ° PF
- Free Motion

- Dorsiflexion Assist
- Double Action, Titanium (Extra Charge)
- Camber Axis
- Oklahoma

Posterior Stop

- Fixed Plastic PF Stop
- Adjustable Snap Stop
- Other (Specify) \_\_\_\_\_

### ACCESSORIES

- Growth Extensions
- Quick Release
- Knee Pad
- 4 Buckle
- 5 Buckle

### KNEE JOINTS

Please Select Option

- Single Medial
- Single Lateral
- Medial/Lateral
- Droplocks **Standard**
- Aluminum **Standard**
- Stainless Steel

Free Motion  
 Posterior Offset

- Ring Lock
- Posterior Offset
  - Ball Catch
  - Lever Lock
  - Ratchet Lock

Release System

- Bail
- BLISS
- HD Lever
- Other (Specify) \_\_\_\_\_