



PRECISION MANUFACTURING

Pivot Point Manufacturing  
3161 Putnam Blvd., Pleasant Hill, CA 94523  
916-775-5712

Pectus

**PATIENT**

Last name: \_\_\_\_\_

First Initial: \_\_\_\_\_

Patient ID: \_\_\_\_\_  N  W

Practitioner: \_\_\_\_\_ Please change to *your* name if necessary.

Facility: \_\_\_\_\_

**Straps:**

Hook & Loop  Dacron Back  
Strap Color, White is Standard \_\_\_\_\_

Transfer Pattern, White is Standard \_\_\_\_\_

**PLASTIC • PADDING**

**Plastic:** 3/16 Thickness is **Standard**  
 PolyPro  CoPoly  HDPE  Other Thickness: \_\_\_\_\_


**Lining: 1/8 thickness is standard**  
 Alplast  EVA  Other Thickness: \_\_\_\_\_

**Lining and Pads:**  
 Chest  Back  Other  
 Rib  Rib Side  L  R

**MEASURE FROM STERNAL NOTCH**

Top Trimline \_\_\_\_\_ Axilla Trimline \_\_\_\_\_ Bottom Trimline \_\_\_\_\_

**SPECIAL INSTRUCTIONS**



- A – M/L at Level of Deformity
- B – Clearly Mark Deformity
- C – Circumference at Level of Deformity
- D – Distance from Apex of Deformity to Sternal Notch

**Please be sure to mark spine and inferior angle of the scapula. Plaster or fiberglass casts are acceptable.**

Please note upon arrival, estimated fabrication is between 3 to 5 days.